

# SPORTident ComCards

## ComCard Pro – Fast. Moscompass Model 2

- o SI-Card 10: 128 punches, 60ms response time
- Needle setting time: 0,5-1 sec.
- Stability during running: good
- Magnet: extra high-power
- Scale: b/w, scale unit 2 grad.

## ComCard Pro – Stable. Moscompass Model 3. (Most Popular)

- o SI-Card 10: 128 punches, 60ms response time
- **Needle setting time:** 1-1,5 sec.
- Stability during running: excellent
- Magnet: extra high-power
- Scale: b/w, scale unit 2 grad

#### **ComCard Pro - Rainbow-Stable.** Moscompass Model 8

- o SI-Card 10: 128 punches, 60ms response time
- Needle setting time: 1-1,5 sec.
- Stability during running: excellent
- Magnet: extra high-power
- Scale: color, scale unit 15 grad.

#### **ComCard Pro - Super-Stable.** Moscompass Model 9

- o SI-Card 10: 128 punches, 60ms response time
- Needle setting time: 1,5-2 sec.
- Stability during running: best
- **Magnet:** extra high-power
- Scale: none

### ComCard Up – Stable. Moscompass Model 11

- o SI-Card 8: 30 punches, 115ms response time
- Needle setting time: 1,5-2 sec.
- Stability during running: good
- Magnet: high-power
- Scale: yellow, 2 grad.